



### Key Features

- Five Alarm Functions
- Drawer Structure
- Strong Refrigeration

### More Detail:

[http://www.haiermedical.com/products/preservation/40\\_deep\\_freezer/dw-40l508.html](http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l508.html)

## Specifications

### Technical Data

Type	-40°C Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	45

### Performance

Cooling performance(°C)	-40
Temperature Range(°C)	-20~-40

### Control

Controller	Microprocessor
Display	LED

## Electrical Data

<b>Power Supply(V/Hz)</b>	220-240/50/208-230/60	
<b>Power(W)</b>	285	430
<b>Electrical Current(A)</b>	2.9	2.9

## Dimensions

<b>Capacity(L/Cu.Ft)</b>	490/17.3	
<b>Net/Gross Weight(approx)</b>	164/200(kg) 361.6/440.9(lbs)	
<b>Interior Dimensions(W*D*H)</b>	685*610*1228(mm) 27.0*24.0*48.3(in)	
<b>Exterior Dimensions(W*D*H)</b>	940*845*1860(mm) 37.0*33.3*73.2(in)	
<b>Packing Dimensions(W*D*H)</b>	997*927*2044(mm) 39.3*36.5*80.5(in)	
<b>Container load (20'/40'/40'H)</b>	12/24/24	

## Functions

<b>High/Low Temperature</b>	Y
<b>Remote Alarm</b>	Y
<b>Power Failure</b>	Y
<b>Sensor Error</b>	Y
<b>Low Battery</b>	Y

## Accessories

<b>Caster</b>	Y
<b>Foot</b>	Y
<b>Test Hole</b>	Y
<b>Drawers</b>	10
<b>Temperature Recorder</b>	Optional
<b>USB Interface</b>	Optional

## CO2 Backup System

## LN2 Backup System

## Chart Recorder

## Others

Certificate

CE