



## Key Features

- CFC-Free Refrigerant
- WHO Approved
- Super Performance

## More Detail:

[http://www.haiermedical.com/products/who\\_pqs\\_certificated\\_equipment/Ice\\_lined\\_refrigerator/hbc-200.html](http://www.haiermedical.com/products/who_pqs_certificated_equipment/Ice_lined_refrigerator/hbc-200.html)

## Specifications

### Technical Data

<b>Type</b>	Ice-lined Refrigerator
<b>Cabinet Type</b>	Chest
<b>Ambient Temperature</b>	5-43°C
<b>Cooling Type</b>	Direct cooling
<b>Defrost Mode</b>	Auto
<b>Refrigerant</b>	CFC-Free
<b>Noise Level (dB)</b>	44

### Performance

<b>Temperature Range(°C)</b>	2-8°C
------------------------------	-------

### Control

<b>Controller</b>	Microprocessor
<b>Display</b>	LCD

### Electrical Data

<b>Power Supply(V/Hz)</b>	220-240/50/208-230/60	
<b>Power(W)</b>	194	226
<b>Electrical Current(A)</b>	1.4	1.53
<b>Power Consumption:stable running(KWh/24h)</b>	1.62	1.5
<b>Power Consumption: cool down test(KWh/24h)</b>	1.52	
<b>Holdover Time at 43°C</b>	more than 25hrs	

## Dimensions

<b>Vaccine Storage Capacity(L/Cu.Ft)</b>	94.5/3.3	
<b>Gross Volume(L/Cu.Ft)</b>	198/7.0	
<b>Net/Gross Weight(approx)</b>	108/114(kg) 238.1/251.3(lbs)	
<b>Interior Dimensions(W*D*H)</b>	970*360*610(mm) 38.2*14.2*24.0(in)	
<b>Exterior Dimensions(W*D*H)</b>	1243*630*818(mm) 48.9*24.8*32.2(in)	
<b>Packing Dimensions(W*D*H)</b>	1325*760*880(mm) 52.2*29.9*34.6(in)	
<b>Container load (20'/40'/40'H)</b>	24/54/81	

## Functions

<b>High/Low Temperature</b>	Y
<b>Power Failure</b>	Y
<b>Sensor Error</b>	Y

## Accessories

<b>Foot</b>	Y
<b>Baskets</b>	9

## Others

<b>Certificate</b>	CE/WHO/PQS/ISO9001/EN ISO13458
--------------------	--------------------------------