



Key Features

- CFC-Free Refrigerant
- WHO Approved
- Super Performance

More Detail:

http://www.haiermedical.com/products/who_pqs_certificated_equipment/Ice_lined_refrigerator/hbc-70.html

Specifications

Technical Data

Type	Ice-lined Refrigerator
Cabinet Type	Chest
Ambient Temperature	5-43°C
Cooling Type	Direct cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise Level (dB)	43

Performance

Temperature Range(°C)	2-8°C
------------------------------	-------

Control

Controller	Microprocessor
Display	LCD

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60	
Power(W)	150	177
Electrical Current(A)	1.1	1.4
Power Consumption:stable running(KWh/24h)	0.64	1.1
Power Consumption: cool down test(KWh/24h)	2.4	
Holdover Time at 43°C	more than 25hrs	

Dimensions

Vaccine Storage Capacity(L/Cu.Ft)	45/1.6	
Gross Volume(L/Cu.Ft)	71/2.5	
Net/Gross Weight(approx)	58/64(kg) 127.9/141.1(lbs)	
Interior Dimensions(W*D*H)	400*360*610(mm) 15.7*14.2*24.0(in)	
Exterior Dimensions(W*D*H)	670*630*818(mm) 26.4*24.8*32.2(in)	
Packing Dimensions(W*D*H)	755*760*880(mm) 29.7*29.9*34.6(in)	
Container load (20'/40'/40'H)	42/90/135	

Functions

High/Low Temperature	Y
Power Failure	Y
Sensor Error	Y

Accessories

Foot	Y
Baskets	3

Others

Certificate	CE/WHO/PQS/ISO9001/EN ISO13458
--------------------	--------------------------------