

-40°C Deep Freezer



	Model	DW-40L92	DW-40L188	DW-40L262	DW-40L348	DW-40L508	DW-40L626	DW-40W100	DW-40W255	DW-40W380				
Technical Data	Cabinet Type	Upright	Upright	Upright	Upright, double door	Upright, double door	Upright	Chest	Chest	Chest				
	Climate Class	N	N	N	N	N	N	N	N	N				
	Cooling Type	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling				
	Defrost Mode	Manual	Manual	Manual	Manual	Manual	Manual	Manual	Manual	Manual				
	Refrigerant	CFC-Free	CFC-Free	CFC-Free	CFC-Free	CFC-Free	CFC-Free	CFC-Free	CFC-Free	CFC-Free				
Performance	Noise(dB)	40	43	44	45	45	49	≤45	≤45	≤45				
	Cooling performance(°C)	-40	-40	-40	-40	-40	-40	-40	-40	-40				
Control	Temperature Range(°C)	-10--40	-10--40	-10--40	-20--40	-20--40	-20--40	-20--40	-20--40	-20--40				
	Controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor				
Electrical Data	Display	LED	LED	LED	LED	LED	LED	LED	LED	LED				
	Power Supply(V/Hz)	220-240/50	220-240/50	220-240/50	115/60	220-240/50	220-240/50	220-240/50	115/60	220-240/50	220-240/50	208-230/60		
	Power(W)	200	290	340	340	360	430	430	255	255	410	550	550	
Dimensions	Electrical Current(A)	2.1	2.5	3.2	4.9	2.9	2.9	2.9	8	5.0	5.0	4.2	5.1	5.1
	Capacity(L/Cu.Ft)	92/3.3	188/6.6	262/9.3	348/12.3	490/17.3	626/22.1	100/3.5	255/9.0	380/13.4				
	Net/Gross Weight (approx)	46/51(kg)	63/71(kg)	88/93(kg)	137/145(kg)	170/200(kg)	298/314(kg)	43/46(kg)	70/82(kg)	82/89(kg)				
	Interior Dimension (W*D*H)	101.4/112.4(lbs)	138.9/156.5(lbs)	194.0/205.0(lbs)	302.0/319.7(lbs)	374.8/440.9(lbs)	657.0/692.2(lbs)	94.8/101.4(lbs)	154/180.8(lbs)	180.8/196.2(lbs)				
	Exterior Dimension (W*D*H)	465*440*645(mm)	380*420*1225(mm)	480*465*1430(mm)	535*607*1228(mm)	685*650*1228(mm)	760*630*1310(mm)	500*385*625(mm)	1040*430*605(mm)	1376*457*625(mm)				
	Packing Dimension (W*D*H)	18.3*17.3*25.4(in)	15.0*16.5*48.2(in)	18.9*18.3*56.3(in)	21.1*23.9*48.3(in)	27.0*25.6*48.3(in)	29.9*24.8*51.6(in)	19.7*15.2*24.6(in)	40.9*16.9*23.8(in)	54.2*18.0*24.6(in)				
	Container load (20'/40'/40'H)	597*610*860(mm)	597*610*1475(mm)	657*685*1665(mm)	750*840*1860(mm)	895*840*1860(mm)	1035*900*1980(mm)	680*600*805(mm)	1243*633*838(mm)	1550*633*838(mm)				
Functions	High/Low Temperature	23.5*24.0*33.9(in)	23.5*24.0*58.1(in)	25.9*27.0*65.6(in)	29.5*33.1*73.2(in)	35.2*33.1*73.2(in)	40.7*35.4*78.0(in)	26.8*23.6*31.7(in)	48.9*24.9*33.0(in)	61.0*24.9*33.0(in)				
	Remote Alarm	Y	Y	Y	Y	Y	Y	Y	Y	Y				
	Hot Condenser	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	Power Failure	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	High/Low Voltage	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	Sensor Error	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	Low Battery	Y	Y	Y	Y	Y	Y	Y	Y	Y				
	High Ambient Temp	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
Accessories	Caster	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	Foot	Y	Y	Y	Y	Y	Y	N/A	N/A	N/A				
	Test Hole	Y	Y	Y	Y	Y	Y	Y	Y	Y				
	Drawers/Inner doors	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	USB Interface	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
Others	Temperature Recorder	N/A	N/A	N/A	N/A	N/A	N/A	Y	N/A	N/A				
	Certificate	N/A	N/A	N/A	N/A	N/A	N/A	Optional	Optional	Optional				