



Key Features

- Drawer Structure
- Strong Refrigeration
- Two Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l92.html

Specifications

Technical Data	
Type	-40°C Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	40
Performance	
Cooling performance(°C)	-40
Temperature Range(°C)	-20~-40
Control	
Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	220
Electrical Current(A)	2.1

Dimensions

Capacity(L/Cu.Ft)	92/3.3
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	435*410*635(mm) 17.1*16.1*25.0(in)
Exterior Dimensions(W*D*H)	640*610*860(mm) 25.2*24.0*33.9(in)
Packing Dimensions(W*D*H)	660*700*935(mm) 26.0*27.6*36.8(in)
Container load (20'/40'/40'H)	54/108/108

Functions

High/Low Temperature	Y
Remote Alarm	N/A
Power Failure	N/A
Sensor Error	Y
Low Battery	N/A
High Ambient Temp	N/A

Accessories

Caster	Y
Foot	Y
Test Hole	N/A
Drawers	3
Temperature Recorder	N/A
USB Interface	N/A

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate

CE
