

Cryo Vials (1.8ml), Internal Thread With Silicone Washer Seal, CAP TOP WITH INSERT



Product Code	Volume (ml)	Size (mm)	Bottom Option	Inner Pack	Case quantity	Cap Insert
GEN39400-2	1.8	12.5x49	Self-Standing	100	1000	White

Product Details:

- # Sterile
- # Made from special polypropylene for storing biological material or cells at temperatures as low as -196 degrees
- # An exclusive silicone washer between cap and vial features a positive seal at any temperature
- # Thick wall makes vials almost unbreakable
- # Large white marking area and graduations in 0.2ml increments.