


# MVE Research Dewars

The MVE Research Dewars are developed specifically for easy and safe transport of liquid nitrogen and samples within facilities. These units are vacuum insulated for superior thermal performance and long hold times.



### Features include:

- All stainless steel construction
- Insulated lid
- Wide mouth for easy access
- Convenient carrying handle (except .5L)
- Meets laboratory safety requirements that prohibit glass lined dewars



	RD-6	RD-3	RD-2	RD-1	RD-1W	RD-0.5
<b>Performance</b>						
LN2 Capacity L	6	3	2	1	1	0.5
Handle	yes	yes	yes	yes	yes	no
<b>Unit Dimensions</b>						
Neck Opening in. (mm)	7.3 (185)	7.3 (185)	3.9 (99)	3.3 (84)	3.9 (99)	2.6 (66)
Inner Diameter in. (mm)	7.3 (185)	7.3 (185)	3.9 (99)	3.3 (84)	3.9 (99)	2.6 (66)
Outer Diameter in. (mm)	7.8 (198)	7.8 (198)	4.8 (122)	4.3 (109)	4.8 (122)	3.4 (86)
Usable Internal Height in. (mm)	10.6 (269)	6.3 (160)	7.1 (180)	6.2 (157)	6.1 (155)	7.1 (180)
Overall Height in. (mm)	11.8 (300)	7.5 (191)	12.3 (312)	9.1 (231)	7.0 (178)	8.0 (203)
NER without cover L/hr	0.4	0.2	0.1	0.1	0.2	0.1

90 Day Warranty

Visit [www.chartbiomed.com](http://www.chartbiomed.com) for more information.

Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Copyright © 2014 Chart Industries Ref ML-CRY00015 Rev C